

"The Clinical Advantage"™

Clear-Lead™ Mobile Barriers

Lead acrylic mobile x-ray barriers

Save Time – new sleek design is easy to clean!



Ionizing Radiation Protection

Clear-Lead™ Mobile Barriers offer large, durable and shatter resistant protection wherever it's needed.

www.biodex.com/barriers

BIODEX

www.biodex.com

1-800-224-6339

Int'l 631-924-9000

CLEAR-LEAD™ MOBILE BARRIERS

Save Time – new sleek design is easy to clean!

Introducing Clear-Lead™ – lead acrylic mobile x-ray barriers, featuring the first design change in thirty years. These barriers are for use around any imaging procedures using ionizing radiation, e.g., fluoroscopy, x-ray, CT, mammography.

The greatest advantage of these barriers over others on the market is the time they save behind the scenes. Their sleek design exposes no hardware, nooks or crannies while the opaque section of the shield is covered with plastic laminate – the result, a smooth, easy-to-clean finish.

Handle-grab edges and large casters make them effortless to maneuver. Clear-Lead™ Mobile Barriers offer large, durable and shatter resistant protection wherever it's needed.



042-580 Clear-Lead Mobile Barrier, Regular Window



042-582 Clear-Lead Mobile Barrier, Full Window

SPECIFICATIONS:

042-580 Clear-Lead Mobile Barrier, Regular Window

- Dimensions: 75" h x 30" w (191 x 76 cm)
Window: 24" h x 30" w (61 x 76 cm)
Opaque Panel: 48" h x 30" w (122 x 76 cm)
Shielded Area: 72" h x 30" w (183 x 76 cm)
- Shielding:
Window: 0.5 mm lead equivalency
Opaque Panel: 0.8 mm lead
- Leg Depth: 10.5" (26.7 cm)
- Casters: Four hospital grade; two locking, two non-locking
- Weight: 102 lb (46.5 kg)

042-582 Clear-Lead Mobile Barrier, Full Window

- Dimensions: 75" h x 30" w (191 x 76 cm)
Window: 60" h x 30" w (152 x 76 cm)
Opaque Panel: 12" h x 30" w (31 x 76 cm)
Shielding Area: 72" h x 30" w (183 x 76 cm)
- Shielding:
Window: 0.5 mm lead equivalency
Opaque Panel: 0.8 mm lead
- Leg Depth: 10.5" (26.7 cm)
- Casters: Four hospital grade; two locking, two non-locking
- Weight: 98 lb (44.5 kg)

042-580 Barrier, Clear-Lead, Regular Window
Window size 24" h x 30" w (61 x 76 cm)

042-582 Barrier, Clear-Lead, Full Window
Window size 60" h x 30" w (152 x 76 cm)

www.biodex.com/barriers

BIODEX

Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9338, Email: info@biodex.com, www.biodex.com