

# Specifications

## Clamps

Clamps are made of anodized aluminum with stainless steel parts. They are durable and rust free.

Standard X-ray table: Fits side of table ranging from 1/2" to 2 1/2" in height.  
Compatible w/ GE, Phillips.

Special Procedures table: Intended for table with C-Arm and no Bucky tray.  
A.Fits side of table ranging from 1 1/2" to 2 1/4" in height.  
B.Fits side of table 3/8"  
All manufacturers.

Radiation Therapy table: Fits side of table ranging from 1/8" to 1 1/4" in height.

## Straps

Straps are adjustable to fit most sizes. They are hook and loop, radiolucent, tear resistant, durable, and machine washable. Straps have "Antimicrobial" backing. Replacement straps are available. Do not put in dryer.

## Custom Design Available

PATENT # D672877